

ULTRA-LIGHTWEIGHT, EASY TO USE AND COMFORTABLE

Compartment barrier gently prevents sliding forward

Equal pressure distribution on the seating surface

Auto-adjusting valve



Air valve automatically adapts the cushion to your body shape and weight

Two compartment layout limits sliding

Easy to clean and disinfect

The Vicair Liberty is an easy to use and comfortable (wheelchair) cushion. It is designed for users who require comfort in combination with a low to medium level of support and pressure ulcer prevention. The cushion is ready-for-use, rain proof and with two compartments it provides an even pressure distribution. This makes it an ideal cushion for occasional use in a large variety of situations.

This Vicair cushion is extremely lightweight.

PRODUCT PERFORMANCE

| | |
|----------------------------------|-------|
| Skin protection / Anti decubitus | ▲▲▲▲△ |
| Pressure redistribution | ▲▲▲▲△ |
| Shear reduction | ▲▲▲▲△ |
| Micro climate control | ▲▲▲▲△ |
| Impact & vibration damping | ▲▲▲▲△ |

Positioning

| | |
|---------------------------|-------|
| Lateral stability | ▲▲▲▲△ |
| Frontal stability | ▲▲▲▲△ |
| Sliding prevention | ▲▲▲▲△ |
| Allows manual adjustments | ▲▲▲▲△ |
| Allows easy transfers | ▲▲▲▲△ |

Maintenance & Safety

| | |
|----------------------------|-------|
| Durability | ▲▲▲▲△ |
| Reliability | ▲▲▲▲△ |
| Functionality when misused | ▲▲▲▲▲ |
| Low maintenance level | ▲▲▲▲▲ |
| Level of hygiene | ▲▲▲▲▲ |

Functionality & Application

| | |
|--------------------------|-------|
| Comfort | ▲▲▲▲△ |
| Low cushion weight | ▲▲▲▲▲ |
| Available size range | ▲▲▲▲△ |
| Maximum user weight | ▲▲▲▲▲ |
| Low seat to floor height | ▲▲▲▲▲ |
| Auto adjustments | ▲▲▲▲△ |

USED WITH

Osteoarthritis/arthritis, rheumatism, spondylitis, pain relief, mobility problems, in a geriatric environment.

TECHNICAL SPECIFICATIONS

| | |
|------------------|-----------|
| Height | 5 cm |
| Average weight | 350 grams |
| Max. user weight | 150 kg |

CUSHION SIZES

- Standard
- Not available

