

Q300 M MINI

Compact in size, BIG on performance.





GO WHERE OTHER POWERCHAIRS CAN'T WITH THE Q300 M MINI.

Ultra-slim at just 520 mm, the all-new Q300 M Mini is the narrowest TRUE mid-wheel drive powerchair - and with a turning circle of just 1200 mm. This chair will take you places that other chairs can't fit.

Narrow doorways?

No problem!

Tight fit in the cafe or elevator?

No Problem.

This 52 cm. base will sneak right in.

The outdoor performance has been taken to the next level too with features like the patented all-wheel suspension - where every wheel is independently cushioned for a safe and comfortable ride. The strong 2-pole motors will power-through the tricky terrain where others would turn back. And there's no need to race back either, the energy efficiency of the Q300 M Mini squeezes up to 30 km. from 34Ah batteries.

So accessing all areas is a breeze, how about accessorising them too? This chair is so good looking, you're always assured a ride in style. We've even upgraded the seating to first-class with the multi-adjustable SEDEO LITE. And with plenty of powered options to choose from, the Q300 M Mini will also cater for your ongoing postural and seating needs.

Got your attention? Okay, let's get started.



COMPACT SPACE? MEET THE ULTRA-SLIM BASE!



ULTRA-NARROW 520 MM WIDTH

Fed up of being told where you can't go?

The Q300 M Mini's ultra-compact 520 mm base is so ridiculously narrow, it'll go where other powerchairs can't. Whether it's a complicated living room, a cluttered kitchen or a narrow doorway, the Q300 M Mini fits in to your life, not the other way around.

520 MM WIDTH

ULTRA-COMPACT 1120 MM TURNING CIRCLE

Thanks to its TRUE mid-wheel-drive technology, the Q300 M Mini turns on the spot, for an ultra-small turning circle of 1120 mm. In comparison, RWD and FWD powerchairs need to 'swing' around on its front or rear drive wheels, so the turning circle is wider.



1120 mm



WHAT IS TRUE MID-WHEEL-DRIVE TECHNOLOGY?

It's when the powerchair's drive wheels are aligned central to its length. The Q300 M Mini provides the most-compact turning circle and increased traction and stability thanks to its TRUE mid-wheel drive design.

SHORT 1060 MM BASE

At just 1060 mm in length, the Q300 M Mini will conveniently fit in-and-out of the smallest of lifts and tight spaces such as kitchens with ease. It's also light enough at 103 kg to be usable with most floor lifts too.





**GIVES YOU SPACE
WHEN YOU NEED IT.**



LOW SEAT HEIGHT FROM 16"

Fit comfortably under tables or when driving into your wheelchair adapted vehicle with the Q300 M Mini's ultra-low seat-to-floor height, starting from 16" (410 mm) (with fixed tilt) and 17" (430 mm) (with powered lift/tilt).



EASY TRANSFERS WITH FLIP-BACK ARMRESTS

If you struggle with transferring in and out of your powerchair, the armrests on the Q300 M Mini easily flip backwards, giving you additional space. You also have a choice of removable or reclining versions too.



GET TIGHTER TURNS WITH THE CENTRE MOUNT LEGREST

Looking for the tightest turn?

The optional centre-mount legrest eliminates the need for traditional legrests at the front of your powerchair. So not only can you make tighter turns without the worry of clipping objects, but you can get closer to tables too. It's also easily adjustable by leg length and angle.



DOES YOUR CONTROLLER GET IN THE WAY?

For those times when you need to get just that little-bit closer, you can swing the controller sideways (or optionally swing the armrest backwards).

THE OUT-PERFORMER.



TRUE OUTDOOR PERFORMANCE

Experience a true compact powerchair with truly sensational outdoor performance. Narrow enough to be nimble in a crowded city, and with the outdoor presence of a much larger powerchair, it's streets ahead of other compact chairs.

ANTI-PITCH TECHNOLOGY

Traditionally when a powerchair comes to a halt on a slope, the stopping motion of the chair will jolt you forwards. That's a problem we've eliminated with Anti-Pitch technology, so you're kept safe and stable in your seat.

ALL-WHEEL INDEPENDENT SUSPENSION

With QUICKIE's all-wheel independent suspension, all six wheels remain on the ground whilst climbing and transitioning obstacles. With its PATENTED suspension, any shocks and bumps are absorbed by the wheel's suspension, so you're assured of a soft and smooth ride.



THE LEAN, KERB-CLIMBING MACHINE.

Don't you hate it when there's no kerb cut in sight? Thankfully the Q300 M Mini with its 12" (300 mm) drive wheels can handle kerbs of up to 3" (75 mm), so there's no need to take the long way round.



NO NEED TO RUSH BACK WITH 30 KM RANGE

Powered by punchy 34 Ah batteries, the Q300 M Mini takes you everywhere where you want to go with an impressive range of up to 30 km .



EASILY NAVIGATE BETWEEN OBSTACLES

The ultra-slim 520 mm base isn't just useful indoors, you can also effortlessly navigate between tightly parked cars or narrow pathways thanks to the Q300 M Mini's ultra-small 520 mm width.

HOW DO YOU ROLL?

GREATER TRACTION

Thanks to its TRUE mid-wheel drive technology, the position of your seat is directly over the drive wheel, so the additional weight is increasing traction onto the wheel and decreasing sliding.

12" DRIVE WHEELS THE INDOOR NAVIGATOR

Interested in the most ultra-compact width for the tightest manoeuvrability?

The Q300 M Mini comes with 12" (300 mm) drive wheels as standard, providing an ultra-narrow width of 520 mm.

DRIVE BASES

YOUR TRAVELLING COMPANION

Travelling on public transport?

The Q300 M Mini's ultra-small width of 520 mm and 1120 turning circle will help you negotiate even the busiest of journeys.





PERFECT FOR HOISTING

Weighing from a slim and trim 103 kg, the Q300 M Mini can be easily hoisted in and out of your vehicle. A real space saver, optional extras such as the fold-down backrest, swing-away legrests or single post removable armrests provide an even tighter, compact package for when space is at a premium.



The Q300 M Mini is crash-test approved for use with Four-Point tie down restraint systems.



RETAIN TRACTION ON RAMPS

Traction can sometimes be lost on ramps as the front castor wheels force the drive wheel off the ground. With the Q300 M Mini, the front castor arms move independently from the base, so the drive wheel will retain its traction.



WHEELCHAIR ADAPTED VEHICLES

Thanks to its ultra-low 16" (410 mm) seat height, the Q300 M Mini will easily fit under the steering wheel of your vehicle. Ultra-narrow, it won't catch on the tracks or side panels whilst in your vehicle either.

YOUR FIRST-CLASS SEAT TO COMFORT.

SEDEO LITE SEATING SYSTEM

Ride in comfort with the SEDEO LITE seating system. Combining an adjustable and configurable seating frame with the broadest choice of support and positioning options, the SEDEO LITE seating is designed to be individually tailored to your needs. Quick and easy to adjust with just one tool, it's the seating system that's easily adaptable to your changing needs.



INTEGRATED HEADREST DESIGN

The integrated headrest mount is compatible with SEDEO and WHITMYER headrests and was created to address a range of user needs from basic support to multiple-point control.



JAY POSTURE & SKIN PROTECTION OPTIONS

Select from a range of optional Jay cushions and backrests, for better comfort and positioning.



TENSION ADJUSTABLE BACKREST

The standard backrest upholstery provides a soft comfort cover with 4 strap tension adjustment to cater for users with a kyphotic posture.



POSITIONING ELEMENTS & INTEGRATED SEAT RAIL

Different positioning elements are available like swing away lateral thoracic supports, thigh supports and knee supports. All of them are easily swing away and detachable. The integrated seat rail offers many options for mounting the different positioning elements.

LENGTH COMPENSATING, SWING AWAY LEGRESTS

The SEDEO LITE's length compensating legrests are easily adjustable to achieve the best possible fit. They also swing-away for easier transfers to and from your Q300 M Mini. Alternatively, the centre-mount legrest allows you sit with your knees bent at a more 90 degree angle, reducing the overall length of your powerchair for tighter manoeuvrability.



JAY & WHITMYER POSITIONING OPTIONS

Do you have a specific clinical seating requirement?

The SEDEO LITE seating system is compatible with a wide range of WHITMYER headrests and JAY cushions and backs are available with the SEDEO LITE seating systems to meet everyone's clinical needs and ensure we offer the best seating possible.

JAY CUSHIONS

Designed for pelvic stability, skin protection and comfort, JAY cushions and backs are low maintenance, highly customisable and are available in a large variety to meet your needs.

WHITMYER HEADRESTS

Created to address a range of user needs from basic support to multiple-point control. WHITMYER headrests come in a variety of pad shapes and sizes for superior head positioning and comfort.



POWER-UP YOUR Q300-M MINI



POWERED OPTIONS

The Q300 M Mini doesn't just have big-powerchair performance, it has the big-powerchair features too! Whether its powered tilt, lift, recline or elevating legrests, the Q300 M Mini is loaded with the options to meet you exact clinical and lifestyle needs.

RELAX AND RISE WITH LIFT & TILT

Reach new heights with 300 mm of seat lift. Providing safer transfers by giving you a higher position to stand from, seat lift also helps you to reach everyday objects that would otherwise be out of grasp. It'll also help to bring you eye-level with the world around you. Combined with 30 degrees of tilt, you can additionally enjoy greater relaxation and pressure distribution on a whole new level.



ENJOY A COMFORTABLE SEATED POSITION WITH TILT-IN-SPACE

Adjust your seated position throughout the day with the Q300 M Mini's 30 degrees of powered tilt. Helping with pressure relief by redistributing your weight from the cushion to the backrest, it'll help to maintain a comfortable seated position with reduced muscular effort .



VR2 JOYSTICK

Every Q400 comes with a reliable VR2 control and a simple to control joystick which allows for intuitive operation.

ATTENDANT CONTROL

Allows for a carer to assist with driving or power seat function control.



VR2 Joystick



Attendant Control

DETAILS THAT SHINE.

LIGHTS AND INDICATORS

Light your way home after a long-day out with the lights and indicators package. They'll keep your journey lit and ensure that you're highly visible to others too.

Now that you have seen all of the Q300 M Mini's features, why not try it for yourself?

Arrange a FREE demonstration today.

TECHNICAL DATA

Q300 M MINI			
Max. user weight:	136 kg	Drive Type	TRUE Mid-Wheel Drive
Seat width:	385 - 510 mm (15" - 20")	Max. safe slope:	6 Degrees
Seat depth:	410 - 560 mm (16 - 22)	Max. kerb climbing:	70 mm (with 12" drive wheels)
Seat height:	410 - 490 mm (16 - 19")	Wheel options:	12"
Backrest height:	460 - 510 mm (18 - 20")	Colour choices:	4 shroud colours
Backrest angle:	-4, 0, 4, 8, 12 Degrees	Max. range (ISO 7176-4)	30 km
Total width:	520 mm (20")	Wheelchair weight:	From 103 kg (including 34 Amp)
Total length:	1060 mm (41.5") - 1180 mm (46.5")	Crash tested:	Yes - approved for forward facing transport in a vehicle (ISO 7176-19). Attachment: 4 chair fixing points with 2 double straps rear & 2 single straps front. No anchored belt restraints allowed. The ease of access and manoeuvrability in vehicles depends on the wheelchair size.
Seat tilt:	0 - 12 Degrees (Fixed Tilt) / 0 - 30 Degrees (Powered Tilt)	Transport without user:	Fix the wheelchair safe in the vehicle (e.g. car or plane)
Turning diameter:	1120 mm (44")	Intended use and environment:	For people with limited mobility, indoor and outdoor use. Wheelchair type class B
Max. speed:	10 kph		
Battery size:	34 Ah		
Electronics:	VR2		

Sunrise Medical Australia Pty Ltd
 11 Daniel Street, Wetherill Park
 NSW 2164
 Phone: (02) 9678 6600
 Email: enquiries@sunrisemedical.com.au
 Website: www.SunriseMedical.com.au

@ 2020 Q300 M Mini Brochure AU Rev A

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. Pictures and illustrations may be overseas models and some seating and accessories not available in Australia. If you are visually impaired, this document can be viewed in PDF format at SunriseMedical.com.au.

